

# Vacuum 660

The Vacuum 660 works in combination with any EMS Physio stimulation unit (ie EMS960, EMS970 and EMS860) to provide easy application of electrodes during treatment. It has been designed to sit underneath the stimulation unit for ergonomic and practical purposes.

- Lightweight, ergonomic and contemporary design
- Pulsed or continuous vacuum
- Powerful vacuum with virtually no vibration
- Supplied with complete set of accessories - plug and play
- Independent cup suction capability
- Suction cup design - allowing sponge insert retention
- Two channel outputs



EMS660



PMA3072 with PMA3071



PMA3073 with PMA3071



**Model : Vacuum 660**

## Technical Specification

**Power Input** 18V DC (via 100-240 Vac 50/60 Hz power supply)

**Vacuum** 0.17 to 0.5 Bar

**Channel Selection** A+B or A only

**Classification** Class 1 Type BF (EN 60601-1)

**Size** 227 x 222 x 55 mm (excluding water trap)  
(D x W x H)

**Weight** 1.5 kg

## Ordering Information

### Vacuum 660

Complete unit includes four 60 mm vacuum electrode cups each with sponge inserts, four vacuum leads and one link lead to connect unit to stimulation unit

**EMS660**

## Accessory Ordering Information

Link lead to Primo Interferential 960  
or Primo Multidyne 970

**PMA3070**

One 60 mm grey vacuum cup

**PMA3071**

One red vacuum lead

**PMA3072**

One black vacuum lead

**PMA3073**

Set of 4 sponge inserts for 60 mm  
vacuum cups

**PMA3074**



Tel : 07736768000



Website : [www.almanara\\_iq.com](http://www.almanara_iq.com)

بغداد / باب المعظم / خلف جامع عاذلة خاتون